

pH Electrodes

pH Meter Accessories

For more help selecting the correct electrode for your application please contact jenwayhelp@bibby-scientific.com

✓ Recommended for application

* Suitable for application

	924 001	027 227	924 067	924 002	924 005	924 003	924 007	924 010	924 015	924 030	924 034	924 047	924 070	924 049	924 050	924 052	924 051	924 076	924 077	924 078	924 079	924 080
Agar																	✓					
Alkalines (high)																		✓			✓	✓
Beer	✓				✓	✓			✓			✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Blood products	*				✓	✓			✓	*		*	✓	✓	*	✓	✓	✓	*	✓	✓	✓
Bread/dough				✓				✓			✓						*					
Cement	*	*		*	✓	*	*	*	*		*	*	*	*	*	*	*	✓	*	✓	✓	*
Cosmetics	✓	*		*	✓	✓	*	*	*		*	*	✓	✓	✓	✓	✓	✓	✓	✓	✓	*
Dairy products	*	*		*	✓	✓	✓	✓	*		✓			✓			✓	✓	*	✓	✓	*
Drinking water	*			*					*			*	*	✓					*	*	*	*
Education	✓	*		*	✓	*	*	*	✓	*	*	✓	*	*	✓	✓	✓	✓	✓	✓	✓	✓
Fats/cream				*	*			✓	*		✓			*	*	*	*	*	*	*	*	*
Field use	✓	✓		✓	*			*				✓			✓	✓	✓	*	✓	*	*	
Fish products	*	*		*	*	*	*	✓	*		✓			*	*	*	✓	*	*	*	*	*
High temperature																					✓	✓
Jam/preserves																			✓			*
Lab environmental							*								✓							
Lab flasks														✓							✓	
Low ionic					✓				*					*	✓			✓		✓	✓	
Meat/cheese				*				✓									✓					
Micro samples						✓	*	*	*	*	*	*	*	*	*	*	*					
Paint						✓		✓						✓			✓					✓
Photographic																		✓				✓
Soil	*	✓		*	*	*	✓	*	*		✓			*	*			*	*	*	*	*
Surface																	✓					
Test tubes							*			✓				✓							*	
Tris buffer										✓					✓							*
Viscous samples				*				*			*						✓		*			