

Datasheet

Ag/AgCl Double Junction Reference Electrode Jenway Part Number 924 017



This electrode is designed for use with any mono electrode and is normally used in the laboratory. This is the best choice for the following special applications:

1. Colloidal suspensions
2. Samples containing iodides
3. Samples that are high in solids
4. Viscous solutions
5. Samples containing proteins

Specifications:

| Parameter | Specification |
|-------------------|-------------------|
| Type | Reference |
| Sealed/Refillable | Refillable |
| End Type | Bulb |
| Temperature Range | 0-50°C |
| Reference Type | AgCl |
| Junction Type | Ceramic Frit |
| Body Material | Glass |
| Dimensions | 120x12mm |
| Cable Length | 1 metre (nominal) |
| Termination | 2mm plug |